

# Animal Health Care

## Urine / Milk Test Analyzer

### UriDoctor VET / MilkDoctor R-50S

#### Specification

- Measurement :** Reflectance Photometer
- Size(Dimension) :** 188(W) x 74(D) x 77(H) mm
- Weight :** 430g
- Power :** DC 12V / 3.33 A or AAA battery 1.5V x 8 pcs
- Throughput :** 45(General mode)-120(Quick mode) tests/hour
- Interface :** USB
- Memory :** 2,000 results
- Display :** 4.3 inch touch screen LCD
- Result printing :** via PC S/W or External printer
- Optional item :** External printer, Bluetooth



#### Available Test

- Urine Test Strip : VET 2AC / VET 10 / VET 11C / VET 11AC
- Milk Test Strip : BHB Milk / LDH Milk / MUN

## VET Urine Test Strip

#### Intended Use

- Diagnosis of veterinary disease for dog, cat, horse, cattle, pig, goat, sheep, etc.

#### Benefits

- Simple and easy to get results
- Including Microalbumin / Creatinine Test

#### Detected disease

- Kidney disease / Bacterial infection / Urinary track infection / Liver disease / Hepatitis / Metabolic Disorder



#### \*Protein : Creatinine Ratio

**The urine protein : creatinine ratio** is a simple test that measures how much protein is being lost through the kidneys.

It is used to determine whether a pet may have serious kidney disease.

It measures whether the excretion of protein is greater than expected when compared to the excretion of creatinine.

The test involves measuring the amount of protein and creatinine in the same urine sample and expressing the difference between the two values as a ratio.



#### Ordering Information

Product List	Uro	Glu	Bil	Ket	S.G	Bld	pH	Pro	Nit	Leu	Ma	Cre
VET - 2AC												
VET - 10												
*VET - 11C												
*VET - 11AC												
Kidney Check DOG/CAT												

## Milk Test Strip for cow

### HEALTH MATE® BHB MILK

**Intended Use :** To monitor nutritional status of dairy cows by checking ketosis

**Measurement Item :** BHB( $\beta$ -hydroxybutyrate)

**Reaction time :** 1 min



### HEALTH MATE® LDH MILK

**Intended Use :** To diagnose udder infection(Mastitis)

**Measurement Item :** Lactate Dehydrogenase(LDH)

**Reaction time :** 2 min



### HEALTH MATE® MUN

**Intended Use :** To diagnose protein nutritional status and feeding program

**Measurement Item :** Milk Urea Nitrogen

**Reaction time :** 1 min



## Urine Ketone Test Strip for cow

### HEALTH MATE® COW KETONE TEST STRIPS

#### Intended Use

- For determination of cow's ketones (Acetoacetate) in urine

#### Benefits

HEALTH MATE® COW KETONE TEST STRIPS is a useful diagnostic tool for determining if individual cows or herds are experiencing elevated ketones. It has been relatively popular due to low cost and ease of conducting the test. The sensitivity and specificity with the urine strip have good results.

